

Detailed Contents

Preface to the Second Edition	xxiii
Preface to the First Edition	xxvi
Acknowledgments	xxviii
Health and Safety	xxix
Intellectual Property	xxx

PART I

Bioconjugate Chemistry

1. Functional Targets

1. Modification of Amino acids, Peptides, and Proteins	3
1.1. Protein Structure and Reactivity	4
Amino Acids	4
Nucleophilic Reactions and the pI of Amino Acid Side Chains	13
Secondary, Tertiary, and Quaternary Structure	15
Prosthetic Groups, Cofactors, and Post-Translational Modifications	19
Protecting the Native Conformation and Activity of Proteins	21
Oxidation of Amino Acids in Proteins and Peptides	23
Solvent Accessibility of Functional Targets in Proteins	29
1.2. Protein Crosslinking Methods	32
2. Modification of Sugars, Polysaccharides, and Glycoconjugates	35
2.1. Carbohydrate Structure and Functionality	36
Basic Sugar Structure	37
Sugar Functional Groups	39
Polysaccharide and Glycoconjugate Structure	44
2.2. Carbohydrate and Glycan Conjugation Methods	49
3. Modification of Nucleic Acids and Oligonucleotides	50
3.1. Polynucleotide Structure and Functionality	51
Nucleotide Functional Groups	53
RNA and DNA Structure	62
3.2. Polynucleotide Crosslinking Methods	66
4. Creating Specific Functionalities	66
4.1. Introduction of Sulphydryl Residues (Thiolation)	67
Modification of Amines with 2-Iminothiolane (Traut's Reagent)	67
Modification of Amines with SATA	71
Modification of Amines with SATP	74

Modification of Amines with SPDP	76
Modification of Amines with SMPT	77
Modification of Amines with N-Acetyl Homocysteine Thiolactone	80
Modification of Amines with SAMSA	81
Modification of Aldehydes or Ketones with AMBH	83
Modification of Carboxylates or Phosphates with Cystamine	84
Modification of Proteins with Cystamine	87
Modification of Nucleic Acids and Oligonucleotides with Cystamine	87
Use of Disulfide Reductants	87
Complete Reduction of Disulfides in Protein Molecules Using DTT	90
Use of DTT to Cleave Disulfide-Containing Crosslinking Agents	91
Ellman's Assay for the Determination of Sulfhydryls	100
4.2. Introduction of Carboxylate Groups	101
Modification of Amines with Anhydrides	102
Modification of Sulfhydryls with Iodoacetate	109
Modification of Sulfhydryls with BMPA*	111
Modification of Hydroxyls with Chloroacetic Acid	113
4.3. Introduction of Primary Amine Groups	114
Modification of Carboxylates with Diamines	114
Modification of Sulfhydryls with N-(β -Iodoethyl)trifluoroacetamide [Aminoethyl-8]	118
Modification of Sulfhydryls with Ethylenimine	119
Modification of Sulfhydryls with 2-Bromoethylamine	120
Modification of Sulfhydryls with 2-Aminoethyl-2'-aminoethanethiolsulfonate	121
Modification of Carbohydrates with Diamines	122
Modification of Alkylphosphates with Diamines	124
Modification of Aldehydes with Ammonia or Diamines	124
Introduction of Arylamines on Phenolic Compounds	125
Amine Detection Reagents	127
4.4. Introduction of Aldehyde Residues	129
Periodate Oxidation of Glycols and Carbohydrates	130
Oxidase Modification of Sugar Residues	131
Modification of Amines with NHS-Aldehydes (SFB and SFPA)	132
Modification of Amines with Glutaraldehyde	134
Periodate Oxidation of N-Terminal Serine or Threonine Residues	136
4.5. Introduction of Hydrazine or Hydrazide Functionalities	139
Modification of Aldehydes with Bis-Hydrazide Compounds	140
Modification of Carboxylates with Bis-Hydrazide Compounds	142
Modification of Amines with SANH, SHNH, or SHTH	143
Modification of Alkylphosphates with Bis-Hydrazide Compounds	146
4.6. Introduction of Saccharide or Glycan Groups	147
Modification of Amines with Mono(lactosylamido) mono(succinimidyl)suberate	149
Modification of Amine or Hydrazide Molecules by Carbohydrates and Glycans	150
Labeling Glycans with Fluorescent 2-Aminopyridine, 2-Amino Benzamide, or Anthranilic Acid	153
Synthesis of Glycosylamines for Conjugating Glycans	155
5. Blocking Specific Functional Groups	156
5.1. Blocking Amine Groups	157

Sulfo-NHS Acetate	157
Acetic Anhydride	158
Citraconic Anhydride	159
Maleic Anhydride	159
5.2. Blocking Sulfhydryl Groups	160
N-Ethyl Maleimide	160
Iodoacetate Derivatives	161
Sodium Tetrathionate	161
Methyl Methanethiosulfonate	163
Ellman's Reagent	164
Dipyridyl Disulfide Reagents	165
5.3. Blocking Aldehyde Groups	166
Reductive Amination with Tris or Ethanolamine	166
5.4. Blocking Carboxylate Groups	167
Tris or Ethanolamine plus EDC	167

2. The Chemistry of Reactive Groups

1. Amine Reactions	169
1.1. Isothiocyanates	170
1.2. Isocyanates	170
1.3. Acyl Azides	171
1.4. NHS Esters	171
1.5. Sulfonyl Chlorides	173
1.6. Aldehydes and Glyoxals	173
1.7. Epoxides and Oxiranes	174
1.8. Carbonates	175
1.9. Arylating Agents	175
1.10. Imidoesters	176
1.11. Carbodiimides	176
1.12. Anhydrides	178
1.13. Fluorophenyl Esters	179
1.14. Hydroxymethyl Phosphine Derivatives	180
1.15. Guanidination of Amines	181
2. Thiol Reactions	182
2.1. Haloacetyl and Alkyl Halide Derivatives	182
2.2. Maleimides	183
2.3. Aziridines	184
2.4. Acryloyl Derivatives	184
2.5. Arylating Agents	185
2.6. Thiol-Disulfide Exchange Reagents	185
Pyridyl Disulfides	186
TNB-Thiol	187

Disulfide Reductants	187
2.7. Vinylsulfone Derivatives	188
2.8. Metal-Thiol Dative Bonds	188
2.9. Native Chemical Ligation	191
2.10. Cisplatin Modification of Methionine and Cysteine	192
3. Carboxylate Reactions	192
3.1. Diazoalkanes and Diazoacetyl Compounds	193
3.2. Carbonyldiimidazole	194
3.3. Carbodiimides	195
4. Hydroxyl Reactions	195
4.1. Epoxides and Oxiranes	195
4.2. Carbonyldiimidazole	196
4.3. N,N'-Disuccinimidyl Carbonate or N-Hydroxysuccinimidyl Chloroformate	196
4.4. Oxidation with Periodate	197
4.5. Enzymatic Oxidation	198
4.6. Alkyl Halogens	198
4.7. Isocyanates	199
5. Aldehyde and Ketone Reactions	200
5.1. Hydrazine Derivatives	200
5.2. Schiff Base Formation	200
5.3. Reductive Amination	201
5.4. Mannich Condensation	201
6. Active Hydrogen Reactions	202
6.1. Diazonium Derivatives	202
6.2. Mannich Condensation	203
6.3. Iodination Reactions	203
7. Photo-Chemical Reactions	204
7.1. Aryl Azides and Halogenated Aryl Azides	204
7.2. Benzophenones	205
7.3. Anthraquinones	205
7.4. Certain Diazo Compounds	207
7.5. Diazirine Derivatives	208
7.6. Psoralen Compounds	208
8. Cycloaddition Reactions	210
8.1. Diels-Alder Reaction	210
8.2. Complex Formation with Boronic Acid Derivatives	210
8.3. Click Chemistry: Cu ¹ -promoted Azide-Alkyne [3+2] Cycloaddition	211

PART II**Bioconjugate Reagents****3. Zero-Length Crosslinkers**

1. Carbodiimides	215
1.1. EDC	216
1.2. EDC plus Sulfo-NHS	219
1.3. CMC	223
1.4. DCC	224
1.5. DIC	226
2. Woodward's Reagent K	228
3. N,N' -Carbonyldiimidazole	228
4. Schiff Base Formation and Reductive Amination	231

4. Homobifunctional Crosslinkers

1. Homobifunctional NHS Esters	235
1.1. DSP and DTSSP	238
1.2. DSS and BS ³	241
1.3. DST and Sulfo-DST	243
1.4. BSOCOES and Sulfo-BSOCOES	244
1.5. EGS and Sulfo-EGS	246
1.6. DSG	248
1.7. DSC	249
2. Homobifunctional Imidoesters	250
2.1. DMA	251
2.2. DMP	252
2.3. DMS	253
2.4. DTBP	254
3. Homobifunctional Sulfhydryl-Reactive Crosslinkers	256
3.1. DPDPB	257
3.2. BMH	258
4. Difluorobenzene Derivatives	259
4.1. DFDNB	259
4.2. DFDNPS	260
5. Homobifunctional Photoreactive Crosslinkers	261
5.1. BASED	262
6. Homobifunctional Aldehydes	262
6.1. Formaldehyde	263
6.2. Glutaraldehyde	265

7. <i>Bis</i> -epoxides	268
7.1. 1,4-Butanediol Diglycidyl Ether	269
8. Homobifunctional Hydrazides	269
8.1. Adipic Acid Dihydrazide	270
8.2. Carbohydrazide	271
9. <i>Bis</i> -diazonium Derivatives	271
9.1. <i>o</i> -Tolidine, Diazotized	272
9.2. <i>Bis</i> -diazotized Benzidine	273
10. <i>Bis</i> -alkyl Halides	274
 5. Heterobifunctional Crosslinkers	
1. Amine-Reactive and Sulfhydryl-Reactive Crosslinkers	277
1.1. SPDP, LC-SPDP, and Sulfo-LC-SPDP	278
1.2. SMPT and Sulfo-LC-SMPT	281
1.3. SMCC and Sulfo-SMCC	283
1.4. MBS and Sulfo-MBS	286
1.5. SIAB and Sulfo-SIAB	288
1.6. SMPB and Sulfo-SMPB	291
1.7. GMBS and Sulfo-GMBS	292
1.8. SIAX and SIAXX	293
1.9. SIAC and SIACX	295
1.10. NPJA	296
2. Carbonyl-Reactive and Sulfhydryl-Reactive Crosslinkers	297
2.1. MPBH	298
2.2. M ₂ C ₂ H	299
2.3. PDPH	300
3. Amine-Reactive and Photoreactive Crosslinkers	302
3.1. NHS-ASA, Sulfo-NHS-ASA, and Sulfo-NHS-LC-ASA	305
3.2. SASD	306
3.3. HSAB and Sulfo-HSAB	308
3.4. SANPAH and Sulfo-SANPAH	310
3.5. ANB-NOS	312
3.6. SAND	312
3.7. SADP and Sulfo-SADP	314
3.8. Sulfo-SAPB	316
3.9. SAED	316
3.10. Sulfo-SAMCA	319
3.11. <i>p</i> -Nitrophenyl Diazopyruvate	322
3.12. PNP-DTP	323

4. Sulfhydryl-Reactive and Photoreactive Crosslinkers	324
4.1. ASIB	325
4.2. APDP	326
4.3. Benzophenone-4-iodoacetamide	328
4.4. Benzophenone-4-maleimide	330
5. Carbonyl-Reactive and Photoreactive Crosslinkers	330
5.1. ABH	331
6. Carboxylate-Reactive and Photoreactive Crosslinkers	332
6.1. ASBA	333
7. Arginine-Reactive and Photoreactive Crosslinkers	333
7.1. APG	334
6. Trifunctional Crosslinkers	
1. 4-Azido-2-nitrophenylbiocytin-4-nitrophenyl ester	336
2. Sulfo-SBED	337
3. MTS-ATF-Biotin and MTS-ATF-LC-Biotin	341
4. Hydroxymethyl Phosphine Derivatives	342
7. Dendrimers and Dendrons	
1. Dendrimer Construction	346
2. Conjugation to Dendrimers	353
2.1. Coupling to Amine-Dendrimers	356
Modification of Amine-Dendrimers with Sulfo-NHS-LC-SPDP	356
NHS-PEG-Maleimide Coupling to Amine-Dendrimers	359
Coupling Glycoproteins to Amine-Dendrimers by Reductive Amination	361
Blocking of Amines on PAMAM Dendrimers	363
Preparation of Sugar-Dendrimer Derivatives	366
Conjugation of Carboxylate Organic Molecules to Amine-Dendrimers	371
Epoxy Activation of Amine-Dendrimers	373
Biotinylation of Amine-Dendrimers	376
Fluorescent Labeling of Amine Dendrimers	380
3. Dendrimer-Chelate Derivatives for Imaging Applications	383
4. Dendrimer Derivatives as Surface Modification Agents	385
5. Dendrimer Fluorescent Quantum Dots	389
8. Cleavable Reagent Systems	
1. Cleavage of Disulfides by Reduction	392
2. Periodate-Cleavable Glycols	393
3. Dithionite-Cleavable Bonds	394
4. Hydroxylamine Cleavable Esters	394
5. Base Labile Sulfones	395

9. Fluorescent Probes

1. Fluorescein Derivatives	400
Amine-Reactive Fluorescein Derivatives	401
Sulfhydryl-Reactive Fluorescein Derivatives	406
Aldehyde/Ketone and Cytosine-Reactive Fluorescein Derivatives	412
2. Rhodamine Derivatives	415
Amine-Reactive Rhodamine Derivatives	416
Sulfhydryl-Reactive Rhodamine Derivatives	425
Aldehyde/Ketone and Cytosine-Reactive Rhodamine Derivatives	427
3. Coumarin Derivatives	430
Amine-Reactive Coumarin Derivatives	431
Sulfhydryl-Reactive Coumarin Derivatives	434
Aldehyde- and Ketone-Reactive Coumarin Derivatives	438
4. BODIPY Derivatives	440
Amine-Reactive BODIPY Derivatives	441
Aldehyde/Ketone-Reactive BODIPY Derivatives	444
Sulfhydryl-Reactive BODIPY Derivatives	449
5. Cascade Blue Derivatives	453
Amine-Reactive: Cascade Blue Acetyl Azide	453
Carboxylate-Reactive: Cascade Blue Cadaverine and Cascade Blue Ethylenediamine	455
Aldehyde/Ketone-Reactive: Cascade Blue Hydrazide	456
6. Lucifer Yellow Derivatives	457
Sulfhydryl-Reactive: Lucifer Yellow Iodoacetamide	458
Aldehyde/Ketone-Reactive: Lucifer Yellow CH	459
7. Phycobiliprotein Derivatives	461
8. Cyanine Dye Derivatives	464
Amine-reactive Cyanine Dyes	467
Thiol-reactive Cyanine Dyes	470
Carbonyl-reactive Cyanine Dyes	472
9. Lanthanide Chelates for Time-resolved Fluorescence	474
10. Quantum Dot Nanocrystals	485
Properties of Quantum Dots	485
Conjugation to QDs	493
Conjugation of Proteins to QDs Using EDC	494
Conjugation to QDs Using sulfo-SMCC	496

10. Bifunctional Chelating Agents and Radioimmunoconjugates

1. DTPA	499
2. DOTA, NOTA, and TETA	500
3. DTTA	501
4. DFA	502
5. Use of Thiolation Reagents for Direct Labeling to Sulfhydryl Groups	503
6. FeBABE	505

11. Biotinylation Reagents

1. Amine-Reactive Biotinylation Agents	507
1.1. D-Biotin and Biocytin	508
1.2. NHS-Biotin and Sulfo-NHS-Biotin	510
1.3. NHS-LC-Biotin and Sulfo-NHS-LC-Biotin	512
1.4. NHS-Iminobiotin	515
1.5. Sulfo-NHS-SS-Biotin	517
2. Sulfhydryl-Reactive Biotinylation Agents	520
2.1. Biotin-BMCC	520
2.2. Biotin-HPDP	522
2.3. Iodoacetyl-LC-Biotin	524
3. Carbonyl- or Carboxyl-Reactive Biotinylation Agents	525
3.1. Biotin-Hydrazide and Biotin-LC-Hydrazide	526
3.2. Biocytin Hydrazide	528
3.3. 5-(Biotinamido)pentylamine	529
4. Photoreactive Biotinylation Agents	530
4.1. Photobiotin	531
4.2. Psoralen-PEO ₃ -Biotin	533
5. Active Hydrogen-Reactive: <i>p</i> -Aminobenzoyl Biocytin, Diazotized	534
6. Glycan Biotinylation Reagents	537
6.1. Biotinylated Aminopyridine	538
6.2. Biotinyl-L-3-(2-naphthyl)-alanine hydrazide	541
6.3. Biotin-PEG-Phosphine	543

12. Iodination Reagents

1. Chloramine-T	548
2. Iodobeads	550
3. Iodogen	553
4. Lactoperoxidase-Catalyzed Iodination	555
5. Iodinatable Modification and Crosslinking Agents	556
5.1. Bolton-Hunter Reagent	557
5.2. Iodinatable Bifunctional Crosslinking Agents	560

13. Silane Coupling Agents

1. Silane Reaction Strategies	565
1.1. Aqueous/Organic Solvent Deposition	566
1.2. Aqueous Deposition	566
1.3. Organic Solvent Deposition	567
1.4. Vapor Phase Deposition	567

2. Functional Silane Compounds	568
2.1. 3-Aminopropyltriethoxysilane and 3-Aminopropyltrimethoxysilane	568
2.2. Carboxyethylsilanetriol	573
2.3. N-(Trimethoxysilylpropyl)ethylenediamine triacetic acid	575
2.4. 3-Glycidoxypolytrimethoxysilane and 3-Glycidoxypolytriethoxysilane	577
2.5. Isocyanatopropyltriethoxysilane	579
14. Microparticles and Nanoparticles	
1. Particle Types	582
2. Particle Characteristics and Stability	584
3. Particle Concentration	588
4. Polymeric Microspheres and Nanospheres	588
4.1. Passive Adsorption	590
4.2. Covalent Coupling to Polymeric Particles	594
4.3. Coupling to Carboxylate Particles	595
4.4. Coupling to Amine Particles	599
4.5. Coupling to Amine Particles Using Crosslinking Agents	599
4.6. Glutaraldehyde	601
4.7. SPDP Coupling to Amine Particles	602
4.8. NHS-PEG _n -Maleimide Coupling to Amine Particles	604
4.9. Coupling to Hydroxyl Particles	606
4.10. Coupling to Hydrazide Particles	613
4.11. Coupling to Epoxy Particles	615
4.12. Coupling to Aldehyde Particles	617
5. Silica Particles	618
5.1. Fluorescent Silica Particles	620
5.2. Silane Functionalization of Silica Particles	625
15. Buckyballs, Fullerenes, and Carbon Nanotubes	
1. Buckyballs and Fullerenes	627
1.1. Properties of Fullerenes	627
1.2. Modification of Fullerenes	629
2. Carbon Nanotubes	638
2.1. Nanotube Properties	638
2.2. Nanotube Functionalization	640
2.3. Detergent or Lipid Modification of Carbon Nanotubes	640
2.4. Pyrene Modification of Carbon Nanotubes	644
2.5. Modification of Carbon Nanotubes by Cycloaddition	645

16. Mass Tags and Isotope Tags

- | | |
|------------------|-----|
| 1. ICAT Reagents | 651 |
| 2. ECAT Reagents | 657 |
| 3. Isobaric Tags | 659 |

17. Chemoselective Ligation: Bioorthogonal Reagents

- | | |
|---|-----|
| 1. Diels–Alder Reagent Pairs | 667 |
| 2. Hydrazine–Aldehyde Reagent Pairs | 669 |
| Protocol for Modification of Amine–Oligo with SANH or SFB | 674 |
| Protocol for Modification of Protein or Antibody with SANH or SFB | 675 |
| Conjugation Using the Aldehyde/Hydrazine Reaction | 675 |
| 3. Boronic acid–salicylhydroxamate Reagent Pairs | 676 |
| 4. Click Chemistry: Cu(1)-promoted Azide–Alkyne [3 + 2] Cycloaddition | 680 |
| 5. Staudinger Ligation | 690 |
| 6. Native Chemical Ligation | 697 |
| 6.1. Expressed Protein Ligation and Inteins | 701 |

18. Discrete PEG Reagents

- | | |
|---|-----|
| 1. Homobifunctional PEG Crosslinkers | 711 |
| 1.1. <i>Bis</i> -NHS Ester PEG Compounds | 711 |
| 1.2. <i>Bis</i> -Maleimide–PEG Compounds | 714 |
| BM(PEG) ₂ , BM(PEG) ₃ , and <i>Bis</i> -MAL–dPEG ₃ | 714 |
| 2. Heterobifunctional PEG Reagents | 718 |
| 2.1. Maleimide–PEG _n –NHS Ester Compounds | 718 |
| 2.2. NHS–PEG _n –Azide/Alkyne Compounds for Chemoselective Ligation | 722 |
| 3. Biotinylation Reagents Containing Discrete PEG Linkers | 726 |
| 3.1. NHS–PEG _n –Biotin Compounds | 727 |
| 3.2. NHS–Chromogenic–PEG ₃ –Biotin | 730 |
| 3.3. Maleimide–PEG _n –Biotin Compounds | 732 |
| 3.4. Hydrazide–PEG ₄ –Biotin | 733 |
| 3.5. Biotin–PEG _n –Amine Compounds | 736 |
| 3.6. Biotin–PEG ₃ –Benzophenone | 739 |
| 4. Discrete PEG Modification Reagents | 739 |

PART III**Bioconjugate Applications****19. Preparation of Hapten–Carrier Immunogen Conjugates**

- | | |
|--------------------------------|-----|
| 1. The Basis of Immunity | 745 |
| 2. Types of Immunogen Carriers | 747 |

2.1. Protein Carriers	748
KLH	748
BSA and cBSA	749
Thyroglobulin and OVA	751
Tetanus and Diphtheria Toxoids	753
2.2. Liposome Carriers	753
2.3. Synthetic Carriers	754
3. Carbodiimide-Mediated Hapten-Carrier Conjugation	755
4. NHS Ester-Mediated Hapten-Carrier Conjugation	763
5. NHS Ester-Maleimide Heterobifunctional Crosslinker-Mediated Hapten-Carrier Conjugation	766
6. Active-Hydrogen-Mediated Hapten-Carrier Conjugation	773
6.1. Diazonium Conjugation	773
6.2. Mannich Condensation	776
7. Glutaraldehyde-Mediated Hapten-Carrier Conjugation	779
8. Reductive Amination-Mediated Hapten-Carrier Conjugation	781
20. Antibody Modification and Conjugation	
1. Preparation of Antibody-Enzyme Conjugates	787
1.1. NHS Ester-Maleimide-Mediated Conjugation	788
Activation of Enzymes with NHS Ester-Maleimide Crosslinkers	789
Conjugation with Reduced Antibodies	790
Conjugation with 2-Iminoethanol-Modified Antibodies	793
Conjugation with SATA-Modified Antibodies	795
1.2. Glutaraldehyde-Mediated Conjugation	797
One-Step Glutaraldehyde Protocol	798
Two-Step Glutaraldehyde Protocol	800
1.3. Reductive-Amination-Mediated Conjugation	800
Activation of Enzymes with Sodium Periodate	802
Activation of Antibodies with Sodium Periodate	803
Conjugation of Periodate-Oxidized HRP to Antibodies by Reductive Amination	804
Conjugation of Periodate-Oxidized Antibodies with Amine or Hydrazide Derivatives	805
1.4. Conjugation Using Antibody Fragments	807
Preparation of F(ab') ₂ Fragments Using Pepsin	807
Preparation of Fab Fragments Using Papain	808
1.5. Removal of Unconjugated Enzyme from Antibody-Enzyme Conjugates	812
Immunoaffinity Chromatography	813
Nickel-Chelate Affinity Chromatography	814
2. Preparation of Labeled Antibodies	816
2.1. Fluorescently Labeled Antibodies	817
2.2. Radiolabeled Antibodies	819
2.3. Biotinylated Antibodies	821

21. Immunotoxin Conjugation Techniques

1. Properties and Use of Immunotoxin Conjugates	827
2. Preparation of Immunotoxin Conjugates	829
2.1. Preparation of Immunotoxin Conjugates via Disulfide Exchange Reactions	833
Pyridyl Disulfide Reagents	834
SPDP	834
SMPT	841
3-(2-Pyridyldithio) Propionate	843
Use of Cystamine, Ellman's Reagent, or S-Sulfonates	844
2.2. Preparation of Immunotoxin Conjugates via Amine- and Sulfhydryl-Reactive Heterobifunctional Crosslinkers	847
SIAB	847
SMCC	850
MBS	852
SMPB	854
2.3. Preparation of Immunotoxin Conjugates via Reductive Amination	855
Periodate Oxidation of Glycoproteins Followed by Reductive Conjugation	855
Periodate-Oxidized Dextran as Crosslinking Agent	857

22. Preparation of Liposome Conjugates and Derivatives

1. Properties and Use of Liposomes	858
1.1. Liposome Morphology	858
1.2. Preparation of Liposomes	861
1.3. Chemical Constituents of Liposomes	863
1.4. Functional Groups of Phospholipids	869
2. Derivatization and Activation of Lipid Components	869
2.1. Periodate Oxidation of Liposome Components	870
2.2. Activation of PE Residues with Heterobifunctional Crosslinkers	871
3. Use of Glycolipids and Lectins to Effect Specific Conjugations	877
4. Antigen or Hapten Conjugation to Liposomes	879
5. Preparation of Antibody-Liposome Conjugates	881
6. Preparation of Biotinylated or Avidin-Conjugated Liposomes	883
7. Conjugation of Proteins to Liposomes	885
7.1. Coupling via the NHS Ester of Palmitic Acid	886
7.2. Coupling via Biotinylated PE Lipid Derivatives	888
7.3. Conjugation via Carbodiimide Coupling to PE Lipid Derivatives	888
7.4. Conjugation via Glutaraldehyde Coupling to PE Lipid Derivatives	890
7.5. Conjugation via DMS Crosslinking to PE Lipid Derivatives	892
7.6. Coupling via Periodate Oxidation Followed by Reductive Amination	893
7.7. Conjugation via SPDP-Modified PE Lipid Derivatives	894
7.8. Conjugation via SMPB-Modified PE Lipid Derivatives	895

7.9. Conjugation via SMCC-Modified PE Lipid Derivatives	896
7.10. Conjugation via Iodoacetate-Modified PE Lipid Derivatives	897
23. Avidin-Biotin Systems	
1. The Avidin-Biotin Interaction	900
2. Use of (Strept)avidin-Biotin Interactions in Assay Systems	902
3. Preparation of (Strept)avidin Conjugates	905
3.1. NHS Ester-Maleimide-Mediated Conjugation Protocols	906
3.2. Conjugation using Periodate Oxidation/Reductive Amination	910
3.3. Glutaraldehyde Conjugation Protocol	913
4. Preparation of Fluorescently Labeled (Strept)avidin	914
4.1. Modification with FITC	915
4.2. Modification with Lissamine Rhodamine B Sulfonyl Chloride	916
4.3. Modification with AMCA-NHS	917
4.4. Conjugation with Phycobiliproteins	918
5. Preparation of Hydrazide-Activated (Strept)avidin	919
6. Biotinylation Techniques	920
7. Determination of the Level of Biotinylation	921
24. Preparation of Colloidal Gold-Labeled Proteins	
1. Properties and Use of Gold Conjugates	924
2. Preparation of Mono-Disperse Gold Suspensions for Protein Labeling	928
2.1. Preparation of 2 nm Gold Particle Sols	928
2.2. Preparation of 5 nm Gold Particle Sols	928
2.3. Preparation of 12 nm Gold Particle Sols	929
2.4. Preparation of 30 nm Gold Particle Sols	929
3. Preparation of Protein A-Gold Complexes	930
4. Preparation of Antibody-Gold Complexes	931
5. Preparation of Lectin-Gold Complexes	932
6. Preparation of (Strept)avidin-Gold Complexes	934
25. Modification with Synthetic Polymers	
1. Protein Modification with Activated Polyethylene Glycols	937
1.1. Trichloro-s-triazine Activation and Coupling	938
1.2. NHS Ester and NHS Carbonate Activation and Coupling	940
1.3. Carbodiimide Coupling of Carboxylate-PEG Derivatives	945
1.4. CDI Activation and Coupling	946
1.5. Miscellaneous Coupling Reactions	948
2. Protein Modification with Activated Dextrans	951
2.1. Polyaldehyde Activation and Coupling	952

2.2. Carboxyl, Amine, and Hydrazide Derivatives	954
2.3. Epoxy Activation and Coupling	956
2.4. Sulfhydryl-Reactive Derivatives	960
26. Enzyme Modification and Conjugation	
1. Properties of Common Enzymes	961
1.1. Horseradish Peroxidase	961
1.2. Alkaline Phosphatase	963
1.3. β -Galactosidase	964
1.4. Glucose Oxidase	965
2. Preparation of Activated Enzymes for Conjugation	966
2.1. Glutaraldehyde-Activated Enzymes	966
2.2. Periodate Oxidation Techniques	966
2.3. SMCC-Activated Enzymes	967
2.4. Hydrazide-Activated Enzymes	967
2.5. SPDP-Activated Enzymes	968
3. Preparation of Biotinylated Enzymes	968
27. Nucleic Acid and Oligonucleotide Modification and Conjugation	
1. Enzymatic Labeling of DNA	970
2. Chemical Modification of Nucleic Acids and Oligonucleotides	973
2.1. Diamine or <i>Bis</i> -Hydrazide Modification of DNA	974
Conjugation via Bisulfite Activation of Cytosine	974
Conjugation via Bromine Activation of Thymine, Guanine, and Cytosine	976
Conjugation via Carbodiimide Reaction with 5' Phosphate of DNA (Phosphoramidate Formation)	978
2.2. Sulfhydryl Modification of DNA	980
Cystamine Modification of 5' Phosphate Groups Using EDC	981
SPDP Modification of Amines on Nucleotides	982
SATA Modification of Amines on Nucleotides	984
2.3. Biotin Labeling of DNA	985
Biotin-LC-dUTP	985
Photobiotin Modification of DNA	987
Reaction of NHS-LC-Biotin with Diamine-Modified DNA Probes	987
Biotin-Diazonium Modification of DNA	989
Reaction of Biotin-BMCC with Sulfhydryl-Modified DNA	990
Biotin-Hydrazide Modification of Bisulfite-Activated Cytosine Groups	990
2.4. Enzyme Conjugation to DNA	992
Alkaline Phosphatase Conjugation to Cystamine-Modified DNA Using Amine- and Sulfhydryl-Reactive Heterobifunctional Crosslinkers	993
Alkaline Phosphatase Conjugation to Diamine-Modified DNA Using DSS	994
Enzyme Conjugation to Diamine-Modified DNA Using PDITC	996
Conjugation of SFB-Modified Alkaline Phosphatase to <i>Bis</i> -hydrazide-Modified Oligonucleotides	998

2.5. Fluorescent Labeling of DNA	998
Conjugation of Amine-Reactive Fluorescent Probes to Diamine-Modified DNA	1001
Conjugation of Sulfhydryl-Reactive Fluorescent Probes to Sulfhydryl-Modified DNA	1002
28. Bioconjugation in the Study of Protein Interactions	
1. Homobifunctional Crosslinking Agents	1006
1.1. DSS and BS ³	1007
1.2. Heavy Atom, Deuterated Crosslinking Agents	1008
1.3. Formaldehyde	1010
1.4. Protein Interaction Reporters	1011
2. Use of Photoreactive Crosslinkers to Study Protein Interactions	1016
2.1. Sulfo-SAND, SANPAH, and Sulfo-SANPAH	1016
2.2. Sulfo-SFAD	1018
3. Trifunctional Label Transfer Reagents	1020
3.1. Sulfo-SBED	1021
3.2. MTS-ATF-Biotin and MTS-ATF-LC-Biotin	1028
4. Metal Chelates in the Study of Protein Interactions	1032
4.1. FeBABE for Protein Mapping Studies	1032
4.2. Ru(II)bpy ₃ ²⁺ for Light-Triggered Zero-Length Crosslinking of Protein Interactions	1037
References	1041
Index	1133